




T&D-scan laser system installations




Drexel University, USA




University of Waterloo, Canada




Shanghai Institute of Technical Physics, China



TU Dortmund, Germany




National Institute of Metrology, China



Saint Petersburg University, Russia

T&D-scan laser system

linewidth of CW radiation is 0,01-0,001 nm in spectral range of 275-1100 nm




wavelength of radiation



computer with control software



Optical reconfigurable laser platform



radiation wavelength meter



analogue-digital control unit

Accessories

optics

circulation units

power meter IR-viewer

